

Automation

REPRESENTED BY
MARCO

“The ability to accurately, quickly, and easily show my patient the differences in the visual experience with their current Rx and their new Rx is incredibly important.”

– John Warren, OD
Warren Eye Care
Racine, WI



MARCO TRS-5100 DIGITAL REFRACTION

- Reduce subjectivity of refractions, and transcription errors
- Standardize reproducible results across all operators
- No more ‘better, 1 or 2’, with immediate comparison of old/new Rx
- Cut testing times in half – with more quality time in patient consults
- Transfer the subjective art of examination to the art of professional interpretation



Accuracy

BECAUSE IT'S TIME FOR BETTER

**Join John Warren, OD | VEE Friday 9:45 AM
1 Hour CE: WaveFront Driven Refraction – The Time Is Now!**

VEE • MS6703



Designed and Manufactured by NIDEK • Represented by Marco • 800.874.5274 • marco.com